

Anterior Chamber Maintainer (Lewicky)

23 Gauge x 9/64"

M2658



23G x 9/64" 0.64 x 3.5mm

Anterior Chamber Maintainer (Lewicky)

150mm Silicone Tubing

Annular rings allow cannula to be self-retaining

Hydrodissector (Angled 90°) (Chang)

27 Gauge x 3/4" 0.5 x 18mm

M3904




27G x 3/4" 0.5 x 18mm Hydrodissector (Angled 90°) (Chang)

Tip flattened in vertical plane

Bimanual Handle, Irrigation

23 Gauge x 5/8" 0.64 x 16mm

M5513




23g x 5/8" (0.64 x 16mm) Cannula

Twin 0.4mm side ports (Curved)

Bimanual Handle, Aspiration

23 Gauge x 5/8" 0.64 x 16mm

M5518



23g x 5/8" (0.64 x 16mm) Cannula
0.3mm side port (Curved)

Coaxial I/A Handpieces

16/21 Gauge


M5580



16/21G 1.6 x 0.81mm I/A Handpiece. Straight tip

16/21 Gauge

M5582



16/21G 1.6 x 0.81mm I/A Handpiece. 45° angled tip

16/23 Gauge

M5585



16/23G 1.6 x 0.64mm I/A Handpiece. Binkhorst tip

Crestpoint Management Ltd.

9100 Watson Road
Suite 100
St. Louis, MO 63126 2241



Phone: (314) 849 7773, (800) 756 7773

FAX: (314) 849 7776

E-mail: sales@crestpointmgt.com

All schematic line drawings, photographs and copy in this Brochure are fully protected by copyright. No part of this Brochure may be reproduced in any form without prior written permission. We reserve the right to make changes at any time, without notice, in product specifications and availability. Schematic and photographic illustrations may vary from actual size. Tip drawings are shown enlarged. Descriptive, typographic, or photographic errors are subject to correction. ZCC402 Revised 4/17/2008 ©2008 Crestpoint Management Ltd.